# For Research Use Only. Not for use in diagnostic procedures.



# Anti-mini-AID-tag mAb-HRP-DirecT

**CODE No.** M214-7

**CLONALITY** Monoclonal

CLONE 1E4

**ISOTYPE** Mouse IgG2a κ

**QUANTITY**  $100 \mu L$ 

**SOURCE** Purified IgG from hybridoma supernatant

**IMMUNOGEN** 17 aa sequence of Auxin Inducible Degron internal region (mini-AID-tag).

FORMURATION PBS/Preservative/Stabilizer

**STORAGE** This antibody solution is stable for one year from the date of purchase when stored at 4°C.

#### APPLICATION-CONFIRMED

Western blotting 1:1,000 for chemiluminescence detection system

REFERENCES 1) Nishimura, K. and Kanemaki, M. T., Curr. Protoc. Cell Biol. 64, 20.9.1-20.9.16 (2014)

2) Kubota, T., et al., Mol. Cell 50, 273-280 (2013)

3) Nishimura, K., et al., Nat. Methods 6, 917-922 (2009)


For more information, please visit our web site <a href="http://ruo.mbl.co.jp/">http://ruo.mbl.co.jp/</a>



RELATED PRODUCTS		PM016	Anti-Luciferase pAb (polyclonal)
Antibodies		PM047	Anti-Renilla Luciferase pAb (polyclonal)
M048-3	Anti-GFP mAb (1E4)	M094-3	Anti-β-galactosidase mAb (5A3)
D153-3	Anti-GFP mAb (RQ2)	PM049	Anti-β-galactosidase pAb (polyclonal)
D153-5 D153-6	Anti-GFP mAb-Biotin (RQ2)	M091-3	Anti-MBP (Maltose Binding Protein) mAb (1G12)
	, <b>2</b> /	M013-3	Anti-Thioredoxin (Trx-tag) mAb (2C9)
D153-8	Anti-GFP mAb-Agarose (RQ2)	PM015	Anti-CBD (Chitin Binding Domain) pAb (polyclonal)
598	Anti-GFP pAb (polyclonal)	PM071	Anti-Calmodulin Binding Protein-tag pAb (polyclonal)
598-7	Anti-GFP pAb-HRP-DirecT (polyclonal)	M211-3	Anti-Strep-tag II mAb (4F1)
PM073	Anti-Renilla GFP pAb (polyclonal)	M214-3	Anti-mini-AID-tag mAb (1E4)
M208-3	Anti-RFP mAb Cocktail (1G9, 3G5)		2 , ,
M155-3	Anti-RFP mAb (8D6)	Smart-IP se	eries
M165-3	Anti-RFP mAb (3G5)	3190	Magnetic Rack
M165-8	Anti-RFP mAb-Agarose (3G5)	M180-11	Anti-HA-tag mAb-Magnetic Beads (TANA2)
M204-3	Anti-RFP mAb (1G9)	M132-11	Anti-HA-tag mAb-Magnetic Beads (5D8)
M204-7	Anti-RFP mAb-HRP-DirecT (1G9)	M185-11	Anti-DDDDK-tag mAb-Magnetic Beads (FLA-1)
PM005	Anti-RFP pAb (polyclonal)	M047-11	Anti-Myc-tag mAb-Magnetic Beads (PL14)
PM005-7	Anti-RFP pAb-HRP-DirecT (polyclonal)	D291-11	Anti-His-tag mAb-Magnetic Beads (OGHis)
M180-3	Anti-HA-tag mAb (TANA2) (200 μL)	D153-11	Anti-GFP mAb-Magnetic Beads (RQ2)
M180-6	Anti-HA-tag mAb-Biotin (TANA2)	M165-11	Anti-RFP mAb-Magnetic Beads (3G5)
M180-7	Anti-HA-tag mAb-HRP-DirecT (TANA2)	M198-9	Anti-E-tag mAb-Magnetic beads (21D11)
561	Anti-HA-tag pAb (polyclonal) (0.1 mL)	M167-11	Anti-V5-tag mAb-Magnetic Beads (1H6)
561-7	Anti-HA-tag pAb-HRP-DirecT (polyclonal)	D058-9	Anti-Multi Ubiquitin mAb-Magnetic beads (FK2)
561-8	Anti-HA-tag pAb-Agarose (polyclonal)	M075-11	Mouse IgG1 (isotype control)-Magnetic Beads
M132-3	Anti-HA-tag mAb (5D8)	M076-11	Mouse IgG2a (isotype control)-Magnetic Beads
M185-3L	Anti-DDDDK-tag mAb (FLA-1) (1 mL)	M077-11	Mouse IgG2b (isotype control)-Magnetic Beads
M185-7	Anti-DDDDK-tag mAb-HRP-DirecT (FLA-1)	M077-11 M081-11	Rat IgG2a (isotype control)-Magnetic Beads
PM020	Anti-DDDDK-tag pAb (polyclonal)	M180-10	Anti-HA-tag mAb-Magnetic Agarose (TANA2)
PM020-7	Anti-DDDDK-tag pAb-HRP-DirecT (polyclonal)	M132-10	Anti-HA-tag mAb-Magnetic Agarose (5D8)
PM020-8	Anti-DDDDK-tag pAb-Agarose (polyclonal)	M132-10 M185-10	Anti-DDDDK-tag mAb-Magnetic Agarose (FLA-1)
M192-3	Anti-Myc-tag mAb (My3) (200 μL)	M047-10	Anti-Myc-tag mAb-Magnetic Agarose (PL14)
M192-6	Anti-Myc-tag mAb-Biotin (My3)	D291-10	Anti-His-tag mAb-Magnetic Agarose (OGHis)
M047-3	Anti-Myc-tag mAb (PL14)	D153-10	Anti-GFP mAb-Magnetic Agarose (RQ2)
M047-6	Anti-Myc-tag mAb-Biotin (PL14)	M165-10	Anti-RFP mAb-Magnetic Agarose (3G5)
M047-7	Anti-Myc-tag mAb-HRP-DirecT (PL14)	M167-10	Anti-V5-tag mAb-Magnetic Agarose (1H6)
M047-8	Anti-Myc-tag mAb-Agarose (PL14)	M198-10	Anti-E-tag mAb-Magnetic Agarose (1110)  Anti-E-tag mAb-Magnetic Agarose (21D11)
562	Anti-Myc-tag pAb (polyclonal) (0.1 mL)	W1196-10	Anti-E-tag mAo-Magnetic Agaiose (21D11)
D291-3	Anti-His-tag mAb (OGHis) (200 μL)	Drotoin Dun	ification Kits
D291-6	Anti-His-tag mAb-Biotin (OGHis)	3320	HA-tagged Protein PURIFICATION KIT
D291-7	Anti-His-tag mAb-HRP-DirecT (OGHis)	3342	HA-tagged Protein Magnetic PURIFICATION KIT
D291-8	Anti-His-tag mAb-Agarose (OGHis)	3342	
D291-A48	Anti-His-tag mAb-Alexa Fluor® 488 (OGHis)		HA-tagged Protein PURIFICATION GEL (1 mL)
D291-A59	Anti-His-tag mAb-Alexa Fluor® 594 (OGHis)	3325	DDDDK-tagged Protein PURIFICATION KIT
D291-A64	Anti-His-tag mAb-Alexa Fluor® 647 (OGHis)	3343	DDDDK-tagged Magnetic Protein PURIFICATION KIT
M089-3	Anti-His-tag mAb (6C4)	3328	DDDDK-tagged Protein PURIFICATION GEL
M136-3	Anti-His-tag mAb (2D8)	2225 205	(5 mL gel)
PM032	Anti-His-tag pAb (polyclonal)	3325-205	DDDDK-tag peptide (1 mg x 5)
PM032-8	Anti-His-tag pAb-Agarose (polyclonal)	3305	c-Myc-tagged Protein MILD PURIFICATION KIT
M167-3	Anti-V5-tag mAb (1H6)	3340	c-Myc-tagged Protein Magnetic PURIFICATION KIT
M215-3	Anti-V5-tag mAb (OZA3)	3306	c-Myc-tagged Protein MILD PURIFICATION GEL
PM003	Anti-V5-tag pAb (polyclonal)	2210	(1 mL gel, 1 mg peptide)
PM003-7	Anti-V5-tag pAb-HRP-DirecT (polyclonal)	3310	His-tagged Protein PURIFICATION KIT
PM003-8	Anti-V5-tag pAb-Agarose (polyclonal)	3311	His-tagged Protein PURIFICATION GEL
PM021	Anti-S-tag pAb (polyclonal)	2215	(1 mL gel, 5 mg peptide)
PM070	Anti-E-tag pAb (polyclonal)	3315	V5-tagged Protein PURIFICATION KIT
PM022	Anti-T7-tag pAb (polyclonal)	3341	V5-tagged Magnetic Protein PURIFICATION KIT
563	Anti-VSV-G-tag pAb (polyclonal)	3316	V5-tagged Protein PURIFICATION GEL (1 mL)
M071-3	Anti-GST-tag mAb (3B2)		
M209-3	Anti-GST-tag mAb (GT5)	Other related antibodies and kits are also available.	
PM022	Anti-GST-tag pAb (polyclonal)	Please vis	it our website at <a href="http://ruo.mbl.co.jp/">http://ruo.mbl.co.jp/</a>
M095-3	Anti-Luciferase mAb (2D4)		
1,10,5 5	The Englishment (221)		

#### **SDS-PAGE & Western blotting**

- 1) Mix the sample with Laemmili's sample buffer, then sonicate briefly (up to 10 sec.).
- 2) Boil the samples for 5 min. and centrifuge. Load the sample per lane in a 1-mm-thick SDS-polyacrylamide gel (12.5% acrylamide) for electrophoresis.
- 3) Blot the protein to a polyvinylidene difluoride (PVDF) membrane at 1 mA/cm<sup>2</sup> for 1 hr. in a semi-dry transfer system (Transfer Buffer: 25 mM Tris, 190 mM glycine, 20% MeOH). See the manufacturer's manual for precise transfer procedure.
- 4) To reduce nonspecific binding, soak the membrane in 10% skimmed milk (in PBS, pH 7.2) overnight at 4°C.
- 5) Wash the membrane with PBS-T [0.05% Tween-20 in PBS] (5 min. x 3 times).
- 6) Incubate the membrane with primary antibody diluted with 1% skimmed milk (in PBS, pH 7.2) as suggested in the **APPLICATIONS** for 1 hr. at room temperature. (The concentration of antibody will depend on the conditions.)
- 7) Wash the membrane with PBS-T (5 min. x 3 times).
- 8) Wipe excess buffer on the membrane, then incubate it with appropriate chemiluminescence reagent for 1 min. Remove extra reagent from the membrane by dabbing with paper towel, and seal it in plastic wrap.
- 9) Expose to an X-ray film in a dark room for 1 min. Develop the film as usual. The condition for exposure and development may vary.



### Western blot analysis of AID (full-length)


Sample: His-tagged AID protein (full-length) + HEK293T cell lysate

Lane 1: 50 ng His-tagged AID protein (full-length)

Lane 2: 25 ng His-tagged AID protein (full-length)

Lane 3: 12.5 ng His-tagged AID protein (full-length)

Immunoblotted with Anti-mini-AID-tag mAb-HRP-DirecT (M214-7)



## Western blot analysis of AID (full-length) after immunoprecipitation

Sample: 1 µg of His-tagged AID protein (full-length) + HEK293T cell lysate

Lane 1: IP with 1 µg of Mouse IgG2a (isotype control) (M076-3)

Lane 2: IP with 1 µg of Anti-mini-AID-tag mAb (M214-3)

Immunoblotted with Anti-mini-AID-tag mAb-HRP-DirecT (M214-7)

Samples were kindly provided by Dr. Masato Kanemaki. (Molecular Function Laboratory, National Institute of Genetics)