

MONOCLONAL ANTIBODY

Early Endosome Marker

Anti-EEA1 mAb

Code No.	Clone	Subclass	Quantity	Concentration
M176-3	3C10	Mouse IgG2a κ	100 μ L	1 mg/mL

BACKGROUND: The early endosome is a cellular compartment inside eukaryotic cells. It is receiving endocytosed material and sorting them for vesicular transport to late endosomes and lysosomes or for recycling to the plasma membrane. EEA1 (Early Endosome Antigen1) is a 170 kDa coiled-coil protein, which is required for vesicular transport of proteins through early endosomes. It binds membrane lipids through its FYVE domain.

SOURCE: This antibody was purified from hybridoma (clone 3C10) supernatant using protein A agarose. This hybridoma was established by fusion of mouse myeloma cell P3U1 with C3H mouse lymphocyte immunized with the synthetic peptide corresponding to N-terminal of human EEA1.

FORMULATION: 100 μ g IgG in 100 μ L volume of PBS containing 50% glycerol, pH 7.2. No preservative is contained.

STORAGE: This antibody solution is stable for one year from the date of purchase when stored at -20°C.

REACTIVITY: This antibody reacts with human EEA1 on Western blotting, Immunoprecipitation, and Immunocytochemistry. The reactivity to mouse and rat EEA1 was confirmed by Western blotting.

APPLICATIONS:

Western blotting: 1 μ g/mL for chemiluminescence detection system

Immunoprecipitation: 2 μ g/200 μ L of cell extract from 3×10^6 cells

Immunohistochemistry: Not tested

Immunocytochemistry: 1 μ g/mL

Flow cytometry: Not tested

Detailed procedure is provided in the following **PROTOCOLS.**

INTENDED USE:


For Research Use Only. Not for use in diagnostic procedures.

REFERENCES:

- 1) Gaullier, J. M., *et al.*, *J. Biol. Chem.* **275**, 24595-24600 (2000)
- 2) Mu, F. T., *et al.*, *J. Biol. Chem.* **270**, 13503-13511 (1995)

SPECIES CROSS REACTIVITY:

Species	Human	Mouse	Rat
Cells	HeLa, A549	NIH/3T3, MEF	NRK
Reactivity on WB	+	+	+



Western blot analysis of EEA1 in HeLa (1), A549 (2), NIH/3T3 (3), MEF (4) and NRK (5) using M176-3.


PROTOCOLS:

SDS-PAGE & Western Blotting

- 1) Wash cells (approximately 1×10^7 cells) 3 times with PBS and resuspend them in 1 mL of Laemmli's sample buffer.
- 2) Boil the samples for 3 minutes and centrifuge. Load 10 μ L of the sample per lane in a 1-mm-thick SDS-polyacrylamide gel and carry out electrophoresis.
- 3) Blot the protein to a polyvinylidene difluoride (PVDF) membrane at 1 mA/cm² for 1 hour in a semi-dry transfer system (Transfer Buffer: 25 mM Tris, 190 mM glycine, 20% MeOH). See the manufacturer's manual for precise transfer procedure.
- 4) To reduce nonspecific binding, soak the membrane in 10% skimmed milk (in PBS, pH 7.2) for 1 hour at room temperature, or overnight at 4°C.
- 5) Incubate the membrane for 1 hour at room temperature with primary antibody diluted with PBS (pH 7.2) containing 1% skimmed milk as suggested in the **APPLICATIONS.** (The concentration of antibody will depend on the conditions.)
- 6) Wash the membrane with PBS-T [0.05% Tween-20 in PBS] (5 minutes x 3 times).
- 7) Incubate the membrane with 1:10,000 of Anti-IgG (Mouse) pAb-HRP (MBL; code no. 330) diluted with 1% skimmed milk (in PBS, pH 7.2) for 1 hour at room temperature.

- 8) Wash the membrane with PBS-T (5 minutes x 3 times).
- 9) Wipe excess buffer off the membrane, and incubate membrane with appropriate chemiluminescence reagent for 1 minute.
- 10) Remove extra reagent from the membrane by dabbing with a paper towel, and seal it in plastic wrap.
- 11) Expose the membrane onto an X-ray film in a dark room for 3 minutes. Develop the film under usual settings. The conditions for exposure and development may vary.

(Positive controls for Western blotting; HeLa, A549, NIH/3T3, MEF and NRK)




Immunoprecipitation of EEA1 from HeLa with mouse IgG2a isotype control (1) or M176-3 (2). After immunoprecipitated with the antibody, immunocomplexes were resolved on SDS-PAGE and immunoblotted with anti-EEA1 polyclonal antibody (MBL; Code no. PM062).

Immunoprecipitation

- 1) Wash cells (approximately 1×10^7 cells) 3 times with PBS and resuspend them in 1 mL of cold Lysis buffer [50 mM Tris-HCl (pH 7.5), 150 mM NaCl, 0.05% NP-40] containing protease inhibitors at appropriate concentrations. Incubate it at 4°C with rotating for 30 minutes; thereafter, briefly sonicate the mixture (up to 10 seconds).
- 2) Centrifuge the tube at 12,000 x g for 10 minutes at 4°C and transfer the supernatant to another fresh tube.
- 3) Add primary antibody as suggested in the **APPLICATIONS** into 300 µL of the supernatant. Mix well and incubate with gentle agitation for 60-120 minutes at 4°C. Add 20 µL of 50% protein A agarose beads resuspended in the cold IP buffer [10 mM Tris-HCl (pH 8.0), 500 mM NaCl, 0.1% NP-40]. Mix well and incubate with gentle agitation for 60 minutes at 4°C.
- 4) Wash the beads 3-5 times with the cold Lysis buffer (centrifuge the tube at 2,500 x g for 10 seconds).
- 5) Resuspend the beads in 20 µL of Laemmli's sample buffer, boil for 3-5 minutes, and centrifuge for 5 minutes. Use 20 µL/lane for the SDS-PAGE analysis.
(See **SDS-PAGE & Western blotting**.)

(Positive control for Immunoprecipitation; HeLa)



Immunocytochemical detection of EEA1 in HeLa using M176-3.

Green: anti-EEA1
Blue: DAPI counter stain

Immunocytochemistry

- 1) Culture the cells in the appropriate condition on a glass slide. (for example, spread 1×10^4 cells for one slide, then incubate in a CO₂ incubator for one night.)
- 2) Wash the glass slide 2 times with PBS.
- 3) Fix the cells by immersing the slide in PBS containing 4% paraformaldehyde for 10 minutes at room temperature.
- 4) Wash the glass slide 3 times with PBS.
- 5) Immerse the slide in PBS containing 0.2% Triton X-100 for 10 minutes at room temperature.
- 6) Wash the glass slide 2 times with PBS.
- 7) Add the primary antibody diluted with PBS containing 2% FCS as suggested in the **APPLICATIONS** onto the cells and incubate for 30 minutes at room temperature (Optimization of antibody concentration or incubation condition are recommended if necessary.)
- 8) Wash the glass slide 2 times with PBS.
- 9) Add 200 µL of 1:500 Alexa Fluor® 488 conjugated anti-mouse IgG (Invitrogen; code no. A11001) diluted with PBS onto the cells. Incubate for 30 minutes at room temperature. Keep out light by aluminum foil.
- 10) Wash the glass slide 2 times with PBS.
- 11) Counter stain with DAPI for 5 minutes at room temperature.
- 12) Wash the glass slide 2 times with PBS.
- 13) Wipe excess liquid off the slide but take care not to touch the cells. Never leave the cells to dry.
- 14) Promptly add mounting medium onto the slide, then put a cover slip on it.

(Positive control for Immunocytochemistry; HeLa)

RELATED PRODUCTS:

- M176-A48 Anti-EEA1 mAb-Alexa Fluor® 488 (3C10)
- M176-A59 Anti-EEA1 mAb-Alexa Fluor® 594 (3C10)
- M176-A64 Anti-EEA1 mAb-Alexa Fluor® 647 (3C10)
- PM062 Anti-EEA1 pAb (polyclonal)
- PM059 Anti-KDEL pAb (polyclonal)
- PM060 Anti-Calnexin pAb (polyclonal)
- PM061 Anti-GM130 pAb (polyclonal)
- M076-3 Mouse IgG2a (isotype control) (6H3)